

Premium Safe Chill Refrigerator

Model SC 290

SC 290 features a rare combination of streamlined and space-saving design with large capacity and state-of-the-art technology. A large control panel displays temperature history, which can be downloaded via the integrated USB port. The exclusive cabinet also features LED lighting, fully electronic temperature control and extractable drawers with partitions for rational storage.



Features

- Large control panel
- Extractable drawers
- LED lighting
- Natural refrigerant
- USB data collection
- Ergonomic full-door handle
- Reversible glass door
- Castors for easy placement
- State-of-the-art design

Specifications for SC 290

Gross volume, litres(cu. ft.)	335 (11.8)
Net volume, litres (cu. ft.)	214 (7.6)
Temperature range, C°	+4°C ±2
Temperature range, at 25° AMB	+4°C ±2
Temperature range, at 32° AMB	+4°C ±2
Power consumption per 24h, kWh	2.49
Rated oad (W)	350
Refrigerant	R600a