

# Premium Ultra Low Temperature -45°C Cabinet

## Model VTS 254

VTS 254 features a rare combination of streamlined and space-saving design with large capacity and state-of-the-art technology. A large control panel displays temperature history, which can be downloaded via the integrated USB port. The exclusive cabinet also features door alarm, high/low alarm and 48h battery backup (optional).



### Features

- Large control panel
- Internal doors
- Data logger
- LED lighting
- Natural refrigerant
- USB data collection
- Ergonomic full-door handle

| Specifications for VTS 254    |            | Dimensions for VTS 254 |               |
|-------------------------------|------------|------------------------|---------------|
| Gross volume, litres(cu. ft.) | 253 (8.9)  | Height, mm (inches)    | 198500 (7818) |
| Net volume, litres (cu. ft.)  | 242 (8.5)  | Width, mm (inches)     | 600 (2362)    |
| Temperature range, C°         | -10 to -45 | Depth, mm (inches)     | 59500 (2342)  |
| Refrigerant                   | Nature R   | Gross weight, kg (lbs) | 135 (297)     |
|                               |            | Net weight, kg (lbs)   | 122 (269)     |