

Low Temperature Freezer -60°C

Model VT 307

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals. Nature R refrigerant, cyclopentane insulation combined with recyclable materials make the freezers extremely eco-friendly.



Features

- Single compressor low-energy system
- Prepared for external logging equipment
- Efficient cyclopentane insulation
- Eco-friendly, HCFC-free Nature R refrigerant
- User-friendly and maintenance-free
- Digital XR30CX controller
- Audible and visual alarm
- Castors

Specifications for VT 307		Dimensions for VT 307	
Gross volume, litres(cu. ft.)	296 (10.5)	Height, mm (inches)	850 (33.46)
Net volume, litres (cu. ft.)	284 (10.0)	Width, mm (inches)	1260 (49.61)
Temperature range, C°	-30 to -60	Depth, mm (inches)	600 (23.62)
Power consumption per 24h, kWh	6,20	Gross weight, kg (lbs)	82 (181)
Rated oad (W)	400	Net weight, kg (lbs)	76 (168)
Refrigerant	Nature R		
No. of baskets	3		