

# Low Temperature Freezer -60°C

## Model 547

The VT low temperature freezers creates the possibility to maintain temperatures as low as -60°C. Supreme stability, reliability, user-friendliness and ease of cleaning make these freezers an ideal solution for laboratories and hospitals. Nature R refrigerant, cyclopentane insulation combined with recyclable materials make the freezers extremely eco-friendly.



### Features

- Single compressor low-energy system
- Prepared for external logging equipment
- Efficient cyclopentane insulation
- Eco-friendly, HCFC-free Nature R refrigerant
- User-friendly and maintenance-free
- Digital XR30CX controller
- Audible and visual alarm

Specifications for VT 547		Dimensions for VT 547	
Gross volume, litres(cu. ft.)	495 (17.5)	Height, mm (inches)	850 (33.46)
Net volume, litres (cu. ft.)	476 (16.8)	Width, mm (inches)	1660 (65.35)
Temperature range, C°	-30 to -60	Depth, mm (inches)	758 (29.84)
Power consumption per 24h, kWh	6,00	Gross weight, kg (lbs)	123 (272)
Rated oad (W)	550	Net weight, kg (lbs)	107 (194)
Refrigerant	Nature R		
No. of baskets	3		