

# VACCINE CARRIERS

## PQS CODE: E004/020, MODEL: BK-VC 2.6 -CF



<b>SPECIFICATIONS</b>	
Vaccine storage capacity	2.6 litres
Weight fully loaded	4.5 kg
Weight empty	1.9 kg
External dimensions	26 x 26 x32 cm
Internal dimensions	19 x 19 x 22. 5 cm
Vaccine storage dimensions	11.5 x 11.5 x 19.5 cm
Lid type & fixing	Removable-
External materials	High density Polyethylene (HDPE)
Internal lining materials	High density Polyethylene (HDPE)
Insulation materials	CFC Free Polyurethane Foam
Insulation thickness	40 MM mm.
Number of icepacks required	4
Number of icepacks supplied	4
Icepack type required	0.6 L compliant with PQS E005/IP01.2
<b>PERFORMANCE</b>	
Cold life without openings	43 hours at 43°C
Warm life	07 hours at -20°C
Cool life	18 hours at 43°C
<b>COMMENTS</b>	
Test Report	Test report: PSB Singapore: S08EEC01362/AT/JT Dt: 28th July 2008
Verification Protocol	PQS E004/VC 01-VP.1:2008
Prices based on	INCOTERMS 2009