

VACCINE CARRIERS
PQS CODE: E004/Under Approval,
MODEL: BK-VC 3.4-CF



SPECIFICATIONS	
Vaccine storage capacity	3.4 litres
Weight fully loaded	6.0 kg
Weight empty	2.55 kg
External dimensions	30.5x30.5x34 cm
Internal dimensions	16 x 16 x 120 cm
Vaccine storage dimensions	8.5x8.5x20 CM
Lid type & fixing	Removable-
External materials	High density Polyethylene (HDPE)
Internal lining materials	High density Polyethylene (HDPE)
Insulation materials	CFC Free Polyurethane Foam
Insulation thickness	40 MM mm.
Number of icepacks required	4
Number of icepacks supplied	4
Icepack type required	0.6 L compliant with PQS E005/IP01.2
PERFORMANCE	
Cold life without openings	45 hours at 43°C
COMMENTS	
Verification Protocol	PQS E004/VC01-VP.1:2008
Prices based on	INCOTERMS 2009