

VACCINE CARRIERS

PQS CODE: E004/022,

MODEL: VDC-24-CF



SPECIFICATIONS	
Vaccine storage capacity	0.9 litres
Weight fully loaded	2.2 kg
Weight empty	1.00 kg
External dimensions	22.5x16x25 cm
Internal dimensions	13 x 10 x 16.5 cm
Vaccine storage dimensions	6x10x16.5 CM
Lid type & fixing	Removable-
External materials	High density Polyethylene (HDPE)
Internal lining materials	Polypropylene (PP)
Insulation materials	CFC Free Polyurethane
Insulation thickness	28 mm.
Number of icepacks required	2
Number of icepacks supplied	2
Icepack type required	0.4 L compliant with PQS E005/IP01.2
PERFORMANCE	
Cold life without openings	20- hours at 43°C
Warm life	06 hours at -20°C
Cool life	03 hours at 43°C
COMMENTS	
Test Report	PSB Singapore: 719169972-EEC-10/01 Dt: 1 St April 2010
Verification Protocol	PQS E004/VC 01-VP.1:2008
2009 PRICES (Incoterms)	