

COMBINED ICELINED REFRIGERATOR/WATERPAKS FREEZER TCW2000 AC



Specifications

Climate zone:	Hot	Min rated ambient temp:	+15°C
Refrigerant:	R134a	Energy source:	Electric Mains 220-240V 50/60Hz or Electric Mains 115/60Hz optional
Appliance tested at:	+43°C	Ext dimensions (HxLxD)	90 x 127 x 78 cm
Performance at:	+43°C	Fuel and cycle type:	Electric - compression
Vaccine storage capacity: (Litres)	60		
Gross volume (Litres):	118	Gross volume:	42
Holdover time (hours):	39.4 hours	Waterpack freezing capacity:	10 kg /24h 12 kg
Energy consumption, stable running (kWh/24 hours):	1.87 kWh/24h	Waterpack storage capacity:	
Energy consumption, cool down test (kWh/24 hours):	4.34 kWh/24h	Energy consumption during freezing:	1.95 kWh/h
Comments:	Electricity source for 115V-60Hz is optional		
Accessories:	4 baskets; 24 waterpacks 0.6L; 4 keys; Documentation		
Spare parts (ref, price):	Electronic Thermostat - 292.2009.00 - 89.36 € Drier - 294.1853.02 - 1.90 €; Fan - 296.9759.17 - 28.30 €; Compressor - 296.9701.41 - 130.32 €; Starting capacitor - 296.9706.48 - 11.54 €		
Shipping volume:	1.05 m3	Shipping weight:	125 kg
Price year	2014	Incoterms	EXW
Price	1-39 units: 2361 €	40-99 units: 2242 €	100 or more units: 2131 €
Quality standard:	- ISO 9001:2008 - ISO 14001 - Other -		
Verification report:	A65998	Verification laboratory:	Intertek UK
Current PQS status:	<i>pre-qualified: 14 Dec 2010</i>	Validity until:	<i>May, 2016</i>