

# Colpostar™ 1H Precision Optics

## 12x Single magnification – accurate, agile and affordable

The ColpoStar 1H offers fiber optic halogen lighting that delivers white, cool, clear natural tone tissue illumination for more accurate examination. Optimal single 12x magnification offers a wider field of view and greater depth of field. ColpoStar 1H provides:

- Distortion-free precision stereo optics
- 300mm working distance
- Green filter
- Stable five-point base with center post design, easily maneuvers
- Fine and coarse height adjustments

### ITEM NUMBER DESCRIPTION

906135 ColpoStar 1H

