

TriScope™ with Trulight™

Overhead suspension arm for full range of motion

The TriScope colposcope offers enhanced stability and the full range of motion needed to complete exams and perform related biopsies and surgery. With fixed triple magnification, TriScope offers a direct view optical pathway to detect even subtle tissue abnormalities.

- Easily adjusts to 8x, 13x and 21x magnifications
- Adjustable overhead suspension arm ensures optimal hand/eye coordination
- Trulight LED provides brighter light without changing bulbs
- High-resolution stereoscopic optical viewing

Accessories

Built-in dual optical beamsplitter attaches digital photographic accessories, video and co-observation teaching tube to the TriScope and TriStar.



Digital USB Video Camera and Adapter



Digital SLR Camera and Adapter

ITEM NUMBER

DESCRIPTION

906150T-4	TriScope with Trulight-4 leg Base
900832	Digital Camera and Adapter Lens for Tri Series Colposcopes
900811U	Digital USB Video Camera for Tri Series Colposcopes
900813	Teaching Tube for Tri-Series Colposcopes




TRISCOPE™
with Trulight™

 LISTED