

# TriStar™ with Trulight™

## Compact unit for space and budget conscious

Our TriStar colposcope features a stable five-point base with center post stand to ensure maximum stability with the same high-resolution optical pathway and fixed triple magnification settings available in the TriScope.

- 8x, 13x and 21x magnification provides clear visualization
- Sturdy, compact base saves valuable space
- Trulight LED lighting offers brighter light without changing bulbs
- High-resolution stereoscopic optical viewing

### Accessories

Built-in dual optical beamsplitter attaches digital photographic accessories, video and co-observation teaching tube to the TriScope and TriStar.



Digital USB Video Camera and Adapter



Digital SLR Video Camera and Adapter

### ITEM NUMBER

### DESCRIPTION

906140T	TriStar with Trulight
900832	Digital Camera and Adapter Lens for Tri Series Colposcopes
900811U	Digital USB Video Camera for Tri Series Colposcopes
900813	Teaching Tube for Tri Series Colposcopes



  
**TRISTAR™**  
with TruLight™

 **LISTED**