

# ZoomStar™ with Trulight™

## Same quality optics on a compact mobile stand

ZoomStar offers the same quality optics and features as ZoomScope on an exclusive compact center stand with a spring loaded main shaft for gross height adjustment, with the addition of fine height adjustment.

- Stereoscopic continuous 4.5 – 20x zoom Nikon optics for optimal viewing and image magnification
- Dual zoom knobs on the optical head allow for left or right handed operation
- Convenient fine focus control handle sharpens the image
- USB Video Camera option

ITEM NUMBER	DESCRIPTION
906057T	ZoomStar with Trulight
906057-40TU	ZoomStar with Trulight Video

