

# WA1000B™ Cryo Console System

**High performance features and professional appearance make the WA1000B™ your best choice.**

The WA1000B™ console offers the ideal combination of easy-to-use controls and built in safety features. It uses the powerful LL100™ two-trigger freezer for Nitrous Oxide. (Carbon Dioxide units are designed differently due to the properties of the gas.) Freeze and defrost only occur when triggered and can be stopped instantly. A console-mounted switch must be activated to permit gas to flow into the LL100™ and a color-coded gauge provides gas pressure reading at a glance. The Temperature Indicator gives surgeons unparalleled control and efficiency during cryosurgery. A large selection of reusable sterilizable cryotips and disposable plastic shields are available.

- All the advantages of the LL100™ Freezer.
- Precise read-out of tip temperature on thermocoupled tips.
- Audible freeze timer starts when freezing temperature is reached and displays total time the tip is actually freezing.

ITEM NUMBER	DESCRIPTION
900506-1-N <sub>2</sub> O	WA1000B Cryosurgical Console – Blank N <sub>2</sub> O LL100
900506-1-CO <sub>2</sub>	WA1000B Cryosurgical Console – Blank CO <sub>2</sub> LL100
900506-2-N <sub>2</sub> O	WA1000B Cryosurgical Console – Digital Temp. (N <sub>2</sub> O) Indicator, Times, Base and LL100
900506-2-CO <sub>2</sub>	WA1000B Cryosurgical Console – Digital Temp. (CO <sub>2</sub> ) Indicator, Times, Base and LLCO <sub>2</sub>

