

Biopsy Punches

For the physician who demands the finest surgical and biopsy instruments, the smart choice is Wallach. We provide a complete collection of the most popular instruments, all handcrafted by seasoned artisans in Europe for exact operating and examining procedures.



Large rectangular bite with upper and lower teeth for maximum sample.

ITEM NUMBER	DESCRIPTION
907028	Tischler-Kevorkian (L 9.3 mm x W 3 mm) (20 cm)
907056	Tischler-Kevorkian Punch (25 cm)



Narrow jaw with additional teeth on the bottom jaw to improve grasping ability and prevent slippage.

ITEM NUMBER	DESCRIPTION
907029	Kevorkian-Pacific (L 8 mm x W 3 mm) (20 cm)
907057	Kevorkian-Pacific Punch (25 cm)



Large bite to assure adequate sample for histological examination.

ITEM NUMBER	DESCRIPTION
907030	Coppleson (L 7.5 mm x W 3 mm) (20 cm)



Small bite size to facilitate sampling of cervical os. Ideal for stenotic cervix.

ITEM NUMBER	DESCRIPTION
907033	Sullivan (L 4 mm x W 2 mm) (20 cm)



Serrated tooth on bottom jaw coupled with unique bite shape to offer precision and control.

ITEM NUMBER	DESCRIPTION
907034	Wittner Straight Standard (L 7.4 mm x W 3.9 mm) (20 cm)

Large bite size and angled jaw improve access to difficult-to-reach areas of the cervix and provide adequate tissue sample.

ITEM NUMBER	DESCRIPTION
907035	Wittner Angled Standard (L 4.4 mm x W 3.9 mm) (20 cm)



Single anchoring tooth and V-shaped jaw facilitate sampling for analysis.

ITEM NUMBER	DESCRIPTION
907036	Appex Standard (L 7 mm x W 3 mm) (20 cm)

Oval cutting jaw provides large tissue sample for analysis.

ITEM NUMBER	DESCRIPTION
907037	Eppendorfer Standard (L 7.4 mm x W 3.9 mm) (20 cm)



Standard Tischler design for minimal discomfort and trauma. Spring handle handcrafted with precision in German stainless steel.

ITEM NUMBER	DESCRIPTION
907008	Tischler Biopsy Forceps (L 7mm x W 3 mm) (20 cm)
907054	Tischler Punch (L 7mm x W 3 mm) (25 cm)



Popular rectangular small tip, Tischler biopsy forceps offers smaller samples with minimal discomfort.

ITEM NUMBER	DESCRIPTION
907014	Baby Tischler Biopsy Forceps (L 4.2 mm x W 2.3 mm) (20 cm)
907058	Mini Tischler (L 4.2 mm x W 2.3 mm) (25 cm)



Teeth are incorporated in a unique Burke-style sharp jaw.

ITEM NUMBER	DESCRIPTION
907031	Burke (L 5.5mm x W 3.25 mm) (20 cm)



Two bottom teeth and one on the top maximizes the ability to grasp the cervix. Larger than Burke-style punch.

ITEM NUMBER	DESCRIPTION
907032	Baggish (L 6.5mm x W 4 mm) (20 cm)

