

Cervical Dilators and Sounds

Craftsmanship and dependability built into every instrument

Small Cervical Dilator Set

This small cervical dilator set is designed for patients with a stenotic cervical os. Set includes 6 dilators ranging from 1 mm to 6 mm in size.

ITEM NUMBER	DESCRIPTION
907049	Small Dilators – Set of 6



Hegar Dilator Set

Stainless steel Hegar dilators provide the option of double-ended dilators in even-millimeter increments, or from 3 to 18 for a greater diameter range over the entire set. German stainless steel. Set of 8 includes: 3/4 mm, 5/6 mm, 7/8 mm, 9/10 mm, 11/12 mm, 13/14 mm, 15/16 mm, 17/18 mm.

ITEM NUMBER	DESCRIPTION
907051	Hegar Dilators – Set of 8 (3/4 mm – 17/18 mm)



Mini Cervical Dilator Set

Ideally suited for nulliparous patients or those with a stenotic os. Includes a universal handle that holds dilator tips from sizes 1 mm to 3 mm in half-millimeter increments. Rounded dilator tips help prevent uterine wall puncture.

ITEM NUMBER	DESCRIPTION
907050	Mini Dilator Set 1 - 3 mm



Sims Uterine Sound

The Sims Uterine Sound is fully graduated for accurate fundal depth determination and features a bulbous tip to reduce the risk of perforation. The Sims Sound is malleable and can adapt to virtually any anatomical variation.

ITEM NUMBER	DESCRIPTION
907048	Sims Sound Grad 12 in



Hank Cervical Dilators

Hank cervical dilators come in a set of 6, ranging from 4.5 mm to 10 mm. Dilators have a full millimeter size on one end and a half millimeter size on the other for smooth dilation.

ITEM NUMBER	DESCRIPTION
907052	Hank Dilators - Set of 6

