

# Rotating Biopsy Punches

**360° Rotating shaft permits biopsies at any angle without obstructing the biopsy site**



Two years of FREE sharpening!  
Standard: 20cm working length  
Long: 25cm working length

Minimal trauma and maximum sample retention are provided by this most-popular bite design.

## ITEM NUMBER

## DESCRIPTION

907004

Tischler Rotating Biopsy Forceps (L 7 mm x W 3 mm) (20 cm)



---

Angled tip facilitates sampling by permitting easier access to difficult areas. Enhances grasping capability.

## ITEM NUMBER

## DESCRIPTION

907019

Angled Tischler Rotating Biopsy Forceps (L 7 mm x W 3 mm) (20 cm)



---

Baby Tischler tip is small enough to permit sampling at the endocervical canal while providing excellent control.

<b>ITEM NUMBER</b>	<b>DESCRIPTION</b>
907059	Rotating Baby Tischler (L 4.2 mm x W 2.3 mm) (25 cm)



---

Smaller size enables sampling of the cervical os in some patients. Permits access to difficult areas.

<b>ITEM NUMBER</b>	<b>DESCRIPTION</b>
907042	Mini-Up Tip (L 4.2 mm x W 2.3 mm) (20 cm)



---

Facilitates sampling of the cervical os and the outer periphery of the ectocervix.

<b>ITEM NUMBER</b>	<b>DESCRIPTION</b>
907039	Mini-Down (L 4.2 mm x W 2.3 mm) (20 cm)



---

Offers improved grasping ability. Prevents slippage of the cervix.

<b>ITEM NUMBER</b>	<b>DESCRIPTION</b>
907040	Kevorkian (L 8 mm x W 3 mm) (20 cm)



---

Moderate bite for reduced trauma.

<b>ITEM NUMBER</b>	<b>DESCRIPTION</b>
907041	Oval (L 5 mm x W 3 mm) (20 cm)



---

An ample, well-oriented sample is provided by the large rectangular bite.

<b>ITEM NUMBER</b>	<b>DESCRIPTION</b>
907060	Kevorkian (L 8 mm x W 3 mm) (25 cm)



---

Reduces tissue trauma. Smaller size enables sampling of the cervical os in some patients.

**ITEM NUMBER**

**DESCRIPTION**

907038

Mini-Townsend (L 4.2 mm x W 2.3 mm) (20 cm)

