

Uncoated Speculum

High quality German engineering from Wallach, a name you can trust

Graves Specula

Traditional design with slightly wider blades.

ITEM NUMBER	DESCRIPTION
903001	Graves Speculum – Small
903002	Graves Speculum – Medium
903003	Graves Speculum – Large
903004	Graves Speculum – Ex-Large



Moore-Graves Speculum

Unique speculum featuring short blades with the same width as a medium Graves.

ITEM NUMBER	DESCRIPTION
903005	Moore-Graves Speculum

Pederson Specula

Thinner blade design.

ITEM NUMBER	DESCRIPTION
903006	Pederson Speculum – Ex-Small
903007	Pederson Speculum – Ex-Narrow
903008	Pederson Speculum – Small
903009	Pederson Speculum – Medium
903010	Pederson Speculum – Large
903011	Pederson Speculum – Ex-Large



Weisman-Graves Specula

Open-sided speculum with one-sided hinge offers easier access for instruments.

ITEM NUMBER	DESCRIPTION
903012	Weisman-Graves Speculum – Right Medium
903013	Weisman-Graves Speculum – Left Medium



View-More Specula

Provides more viewing room at the yoke and allows for easier access for instruments.

ITEM NUMBER	DESCRIPTION
903014	View-More Speculum – Medium
903015	View-More Speculum – Large



View-Maxi Speculum

Extra-wide rectangular yoke for enhanced viewing and easier access.

ITEM NUMBER	DESCRIPTION
--------------------	--------------------

ITEM NUMBER	DESCRIPTION
903016	View-Maxi Graves Speculum – Medium



Omni Specula

Offers narrower blades combined with lateral blades that lock in place for better visualization.

ITEM NUMBER	DESCRIPTION
903028	Omni Speculum – Medium 4 Blade
903029	Omni Speculum – Large 4 Blade

