

Uncoated Wall Retractor and Uncoated Endocervical Speculum

Uncoated precision German stainless steel instruments

Uncoated Wall Retractor

Contoured and shaped to fit into the yoke of a standard speculum, the retractor expands beyond the width of the speculum blades and retracts the side walls of the vaginal vault to achieve an unobstructed view of the cervical anatomy.

ITEM NUMBER	DESCRIPTION
907043	Lateral Retractor View-More Opening



Uncoated Endocervical Speculum

Standard ratchet handle. Choose from a standard tip for use on multiparous women or a narrow tip particularly useful in nulliparous or small women, post menopausal patients, or those with a stenotic os.

ITEM NUMBER	DESCRIPTION
907010	Kogan Endocervical Speculum (5mm tip)
907015	Kogan Endocervical Speculum – Narrow (3mm tip)

