

Electrosurgical System Accessories

Accessories for Wallach's Quantum 2000™ and other comparable products

Hand Switch Operated Pencils

For use with Quantum 2000 and other compatible electrosurgical generators. Sterile and single use only. Convenient cut and coag push buttons.


ITEM NUMBER	DESCRIPTION
909089	Hand Switch Operated Pencils, Sterile (25 per Box)



Disposable Patient Split Grounding Pads (with Connector)

For use with Quantum 2000, these pads offer extra patient safety by signaling the system to shut down if the patient is not in complete contact with the pad.


ITEM NUMBER	DESCRIPTION
909078	Disposable Patient Split Grounding Pads with Connector (10 per Box)



Foot Pedal Operated Pencils

For use with Quantum 2000 and other comparable electrosurgical generators. Sterile and single use only. Reusable adapter required. Use with any standard disposable electrodes.

ITEM NUMBER	DESCRIPTION
909028	Foot Pedal Operated Pencils, Sterile (25 per Box)



Ferric Sub sulfate Aqueous Solution

Make postoperative bleeding easier to control. Twelve single-dose 15 ml bottles come in a convenient box.

ITEM NUMBER	DESCRIPTION
909084	Ferric Sub sulfate Aqueous Solution (12 per Box)



Electrosurgical Disposable Grounding Pads (for use with reusable cable)

Universal, cost-effective electrosurgical split grounding pad for use with Quantum 2000 and other comparable units. Must be used with durable, reusable 9-ft 4-in cable that can be disinfected after each use.

ITEM NUMBER	DESCRIPTION
909152	Electrosurgical Disposable Grounding Pads for Use with Reusable Cable (5 per Box)
909153	Reusable Cable



Dental Style Anesthetic Cartridge Syringe

Reusable syringe accepts lidocaine cartridges for quick and easy anesthetic blocks.

ITEM NUMBER	DESCRIPTION
909171	Dental Style Cartridge Syringe



Disposable Smoke Evacuation Speculum Tubing (Non-Sterile)

Speculum Tubing (Non-Sterile) Slides through slot of compatible specula and attaches to flexible smoke evacuation tubing.

ITEM NUMBER	DESCRIPTION
920001	Disposable Speculum Tubing (50 per box)



Endocervical Block Kit (for Intra-Cervical Block)

- Single use only, reinforced, long EndoBlock needle
- 18 gauge needle for drawing anesthetic
- Disposable syringe with finger grip

ITEM NUMBER	DESCRIPTION
-------------	-------------

ITEM NUMBER	DESCRIPTION
920009	Endocervical Block Kit (10 Kits per Box)



QKIT Procedure Tray for Loop Electrosurgical Excision Procedure

- Two cartridges containing 2% Lidocaine with Epinephrine
- Endocervical block needle, long threaded hub
- Lugol's (strong iodine solution)
- Ferric subsulfate aqueous solution
- Underpad, panty shield, alcohol wipe, and six large cotton swabs

ITEM NUMBER	DESCRIPTION
909074	QKIT Procedure Tray (5 Trays per Box)



EndoBlock Needle (for Intra-Cervical Block)

- Sterile, single use only, reinforced and extended shaft
- Threaded hub needle for dental style cartridge syringe
- Luer Lock for use with medical syringe

ITEM NUMBER	DESCRIPTION
920017	EndoBlock Needle Luer Lock (6 per Box)
920017-1	EndoBlock Needle Luer Lock (60 per Box)
920021	EndoBlock Needle Threaded Hub (6 per Box)



Smoke Evacuator Filter Assembly

Filter and evacuation tubing efficiently removes airborne vapor that may contain potentially harmful microbes. Maximum operating time of 30 minutes.

ITEM NUMBER	DESCRIPTION
909094	Filter Assembly for Biovac Smoke Evacuator (5 per Box)



Disposable Smoke Evacuation Tubing with Attached Speculum Tubing

Six-foot long sterile, single use only smoke evacuation tubing with attached speculum tubing connects to the smoke evacuator filter assembly.

ITEM NUMBER

DESCRIPTION

920002

Disposable Smoke Evacuation Tubing with Speculum Tubing (10 per box)

