

Disposable Loop Electrodes

Disposable, extra-rigid loop and ball electrodes meet a wide variety of needs

- Sterile, single use only
- Tungsten wire
- 13cm shaft (unless otherwise specified)
- Five per box
- High quality, made in the U.S.A.



3mm Ball

909038



5mm Ball

909003



5mm Ball

(6cm Shaft)
909136



10mm x 7mm Loop

(1.0cm x 0.7cm)
909133



15mm x 5mm Loop

(1.5cm x 0.5cm)
909134



15mm x 8mm Loop

(1.5cm x 0.8cm)
909011



10mm x 10mm Loop

(1.0cm x 1.0cm)
909007



20mm x 8mm Loop

(2.0cm x 0.8cm)
909013



0.8mm Dia. X 16mm Needle

(5.5cm Shaft)
909135



10mm x 10mm Square

(1.0cm x 1.0cm)
909131



0.8mm Dia. X 25.4mm Large Needle

(12.7cm Shaft)
909005



10mm x 8mm Square

(1.0cm x 0.8cm)
909017



20mm x 10mm Loop

(2.0cm x 1.0cm)
909132



20mm x 15mm Loop

(2.0cm x 1.5cm)
909009