

Medical Emergency Cart

Model CA3000-SL



Product features

- Mild steel plate construction
- Stainless steel work top
- Scratch & chip resistant epoxy powder powder coating
- Collapsible work tray
- 4 drawers with sealed lock and ABS molded bins in top 2 drawers
- Cardiac board located at back of cart
- Equipped with protective bumper guard rail
- Provision for Oxygen tank & IV pole
- Mobile on 4 x 100mm dia. swivel castors with 2 brakes

Technical Specification

- Dim (mm): 660(W) X 460(D) X 890(H)
- Optional : Central locking with key (CA3000)