

Medication Cart

Model CA1150



Product features

- Aluminium-alloy, ABS and stainless steel construction
- Consist of 6 layers of drawers designed for dual-sided usage with self-locked mechanism when closed.
- Mobile on 4 x 125mm dia. swivel castors with 2 brakes
- Provide with 50 medicines boxes
- Accessories inclusive of:- - 2 x Litter buckets - 1 x Syringe disposal - 1 x Adjustable board - 1 x Transparent file box 1 x Set of 5-hooked hanging device

Technical Specifications

- Dim(mm): 850(W) x 520(D) x 1010(H)
- Optional Accessories: Dispenser & Shelve