

Examination Table for Gynae

Model EC7000(GS)



Product features

- Tubular steel leg frame with heavy base (3 section)
- Scratch & chip resistant epoxy powder coating
- Gas spring backrest and cut-away/lowered-down footrest
- Height adjustable by hydraulic system foot pedal at both sides of bed
- Upholstery with Mediflex cover with PU foam of 75mm thick
- Inclusive of a pair of adjustable legrest , a stainless steel tray, and Paper roll holder at head end
- Mobile on 4 x 125mm swivel castors with 2 brakes

Technical Specification

- Dimension (mm) : 1880(L) x 610(W) x 610-960(H)
- Optional accessories : tilting + telescopic siderail