

Hospital Hydraulic Bed

Model BA300-T



Product features

- Steel bed frame with scratch and chip resistant baked epoxy powder powder coating
- Hilo adjustable by hydraulic jack controlled by foot pedal
- Backrest & kneerest adjustable by gas spring mechanism
- Trendelenburg/reverse trendelenburg assisted by gas spring
- Detachable ABS head and foot boards
- 125mm dia. Castor with 2 brakes diagonally fitted
- Inclusive of collapsible siderail and mattress
- provision of 4 I.V Drip Pole attachment
- Optional : Central locking castor

Technical Specifications

- Safety load capacity: 250kg
- Hilo: 460 - 820mm
- Backrest: 0°- 80°
- Kneerest: 0°- 40°
- Trendelenburg: 0°- 12°
- Reverse Trendelenburg: 0°- 6°
- Overall Dim (mm): 2140(L) x 1040(W)