

# Hospital ICU/CCU Electrical Bed

## Model BA8989



### Product features

- Steel bed frame with scratch and chip resistant baked epoxy powder coating
- PE Head & Foot boards with locking mechanism for easy attachment, detachment and can be cleaned easily
- Split fold down PE safety side rails with one-hand operation
- Backrest with X-Ray radio Translucent base complete with X-Ray cassette tray
- Perforated Mattress support base
- Four high-quality 150mm diameter German Twinwheel with central locking system( free, steer and total brakes)

- Hi-Lo, Backrest, Kneerest, Trendelenburg and Reverse Trendelenburg by electrical motorized
- Inner and outer sides of backrest side rails embedded with control buttons
- Control panel Embedded onto the backrest side rails internally, and externally
- Nurse lock and nurse handheld control located at foot end
- Nurse lock to lock out total or individual function of control buttons
- Nurse handheld control with CPR button to neutralize bed position during emergency
- Quick CPR release lever for flattening the backrest located at both sides of the bed
- Auto contour function to avoid patients slipping down while performing cardiac chair position
- 6 provisions for IV pole attachment
- 2 provisions for monkey pole
- Drainage bag holder provided
- PU bumper at 4 corners
- IV drip pole storage

### **Technical Specifications**

- Backrest: 0° - 70°
- Hi-Lo: 480-780mm
- Legrest: 0° - 35°
- Kneerest: 0° - 20°
- Trendelenburg & Reverse Trendelenburg: -15o to +15o
- Auto contour
- Cardiac Chair Position
- Vascular Position
- Overall Dim(mm): 2245(L) x 1090(W)
- Safety load capacity: 250kg UDL

### **Electrical System**

- Linak DC-actuator system (TUV & UL approval & EC60601-1)
- 24 VDC with built in back-up battery
- Input Voltage: 220 - 240V; 50Hz - 60Hz
- Safety Rating Control box: IPX4,IPX6 Actuators: IPX4, Nurse Control: IPX6

### **Accessories Included**

- Mediflex mattress with PU foam x 1 unit Mattress Dim(mm):1980(L) x 915(W) x 100(H)
- Mattress cover with fire retardant, anti-bacterial, anti-static and fluid resistant x 1 unit
- Telescopic IV Drip Pole with 2 hooks (Stainless steel) x 1 unit

### **Optional Accessories**

- Monkey pole
- Oxygen cylinder holder
- Portable monitor/Defibrillator tray
- Orthopedic Frame
- Patient Mate Mattress with memory foam
- Pressure reducing mattress