

Hospital Fixed Height Bed

Model BA2010-P



Product features

- Steel construction
- Scratch & chip resistant epoxy powder coating
- backrest & kneerest adjustable by manual crank winding mechanism at foot end of bed
- Easy detachable PE head & leg boards with lock
- built in bumper
- provision of IV drip pole attachment
- mobile on 4 x 100mm swivel grey rubber wheels diagonally fitted with brakes *Available in single fowler for backrest adjustable only (BA3010-P)

Technical Specifications

- Backrest: 0° - 80°
- Overall Dim (mm): 2110(L) x 955(W) x 570(H)