

# Child Cot

## Model DC2000-C



### Product features

- Steel base with ventilation holes
- Available with wooden head/head board with cartoon picture
- Available with wooden head/head board with cartoon picture
- Scratch and chip resistant epoxy powder coating
- Safety side rails of filling not exceeding 102mm spacing with rubber bushes along 2 chrome rods fitted on leg frame eliminating friction noises
- Auto engage side rail
- Mobile on 4 x 100mm swivel castors with diagonally fitted brakes

### Technical Specifications

- Dimension (mm): 1400(L) x 850(W) x 665(H) (from floor to mattress platform)
- Mattress platform to top of side rail: 590mm (H) - *\*optional accessories: Mediflex mattress*
- Mattress Dimension (mm): 1285(L) x 762(W) x 102(Thick) - *\*available for adjustable backrest by mechanical crank winding mechanism in model DC2010-C (Single fowler)*