

Patient Transport Trolley

Model PT3300-GS



Product features

- Hydraulic foot pedal controlled Hi-Lo position
- Tilting assisted by gas spring
- Backrest adjustable by gas spring
- PVC bumper guard rail
- Push handle at head and foot ends
- Collapsible side rails
- IV pole
- Storage utility tray
- Oxygen tank holder

Technical Specifications

- Material : Mild steel frame work
- Finishing : Epoxy powder coating
- Safety load : 250kg
- Backrest : 0° - 90°
- Trendelenburg : 0° - 12°
- Reverse Trendelenburg : 0° - 6°
- Castor Dim (mm) : 200mm (central locking)
- Mattress Dim (mm) : 1930(L) x 660(W) x 50(H)
- Approx. Dim (mm) : 2100(L) x 787(W) x 558 - 878(H)

