

Patient Transport Trolley

Model PT2000



Product features

- Removable / interchangeable head panel
- Sliding epoxy side rails
- Provision for I.V pole and oxygen cylinder
- Protective bumper located at four corners
- Storage utility tray
- Trendelenburg / tilting and backrest positioning by means of smooth winding mechanism

Technical Specifications

- Material : Mild steel frame work
- Finishing : Epoxy powder coating
- Castor Dim (mm) : 125mm dia. (*optional: 200mm dia.)
- Approx Dim (mm) : 2083(L) x 711(W) x 794(H)