

1600 Auto-Chemistry Analyzer



- > Holographic concave flat field grating with rear spectrophotometry and cluster optical path to achieve micro volume detection of reaction solution.
- > Large-volume sample disk to provide more positions
- > 60nm polished probe with functions of clot detection, collision protection and automatic resetting.
- > Liquid level detection to reduce carry-over contamination on probe surface.
- > 2 refrigerated reagent compartments with 134 reagent positions.
- > Simultaneous analysis of 3 ISE items.
- > Connectable to LIS/HIS system to realize remote operation and maintenance.