

300B Auto-Chemistry Analyzer



The CS-300B auto-chemistry analyzer is a small-sized floor top analyzer. It is the ideal instrument for consolidating testing in small or medium size laboratories with low workload. It offers a performance of up to 66 photometric items and 3 ISE items analysis simultaneously. And STAT is achieved through very easy operation.

- > Holographic concave flat field grating with rear spectrophotometry and cluster optical path to achieve micro volume detection of reaction solution.
- > Large-volume sample disk to provide more positions
- > 60nm polished probe with functions of clot detection, collision protection and automatic resetting.
- > Liquid level detection to reduce carry-over contamination on probe surface.
- > 67 reagent positions and 115 sample positions.
- > Constant throughput is 300T/H and ISE 300 T/H.
- > Connectable to LIS/HIS system to realize remote operation and maintenance.