

T240 Auto-Chemistry Analyzer



The CS-T240 auto-chemistry analyzer offers the ultimate in convenience, performance and confidence for small to medium size laboratories. This automated random access clinical chemistry analyzer can be easily integrated with any type of laboratories for routine, emergency analysis

- > Holographic concave flat field grating with rear spectrophotometry and cluster optical path to achieve micro volume detection of reaction solution.
- > Large-volume sample disk to provide more positions
- > 60nm polished probe with functions of clot detection, collision protection and automatic resetting.
- > Liquid level detection to reduce carry-over contamination on probe surface.
- > Shared disk of reagent and sample with total 66 positions.
- > Constant speed 240T/H.
- > Connectable to LIS/HIS system to realize remote operation and maintenance.