

Ampla 2085



- Optional items that allow configurations to the ICU, post-surgery or baby nursery;
- Three options for screens: LED display, monochrome and color LCD;
- Three options for bed – electric or manual table and acrylic bassinet;
- Bypass drawer for X-rays;
- Optional for gas panel with aspirator and/or resuscitator unit;
- Two options for coupled phototherapy: Bilitron or Bilitron Bed;
- Two options for warming: radiant warmer – articulate reflector – or thermal mattress;
- Optional kit for transport, which includes thermal mattress and 2-hour batteries;
- Possibility of coupling Bubble CPAP – Babypap;
- In-bed scale and pulse oximeter.