

SUCTION PUMP (HI-VAC JR)



Technical features

Housing	MS powder coated and four castors
Max. vacuum/LPM	700 mm Hg.50ltrs/min
Pump type	Double rotary vane pump
Jars	2x2 ltrs PC jar
Filter	Trap filter
Tubing	10mm ID× 2 mtr.(PVC)
Vacuum gauge	10cm dia 0-760 mm Hg
Power supply	220 v AC, 50/60 hz,380 watt
Noise level	50 db A ± 3
Dimensions	49×36×70 cms
Net weight	37 kgs
Overflow safely	Mechanical
Accessories	110 v on request
Certifications	CE,IEC 60601-1